

BRIEF COMMUNICATION

P SYSTEM ANTIGENIC DETERMINERS EXPRESSION IN *Ascaris lumbricoides*

Patricia PONCE DE LEÓN & Juana VALVERDE

SUMMARY

The P System antigens have been detected in numerous parasites, bacterias and viruses, nevertheless the clinical significance is still unknown. The aim was to study the presence of P₁ antigenic determiners in *A. lumbricoides* extracts by means of the use of 6 different monoclonal antibodies of well-known concentrations and Ig class. We worked with 14 *A. lumbricoides* extracts. Inhibition Agglutination Test was made in a bromelin enzymatic medium and 4 °C temperature. Titre, Score and Sensitivity Parameter were determined for each monoclonal antibody against red cells suspension used as revealing system. Ten extracts inhibited the agglutination of all anti P₁ monoclonal antibodies. The 4 remaining extracts only inhibited the agglutination of some of them. It is demonstrated that the extracts have P₁ activity. This activity is independent of titre, Score, Sensitivity Parameter, concentration and Ig class and it depends on the epitope at which the monoclonal antibody is directed.

KEYWORDS: *Ascaris lumbricoides*; P System

The P blood group System was discovered by LANDSTEINER & LEVINE in 1927⁴ in the course of a systematic attempt to identify new alloantigens by the immunization of rabbits with human erythrocytes.

The new antigen was named P, but it was later designated P₁ when the antigen now known as P was identified^{5,14}. The serology of this system has gradually increased in complexity, and the recognition of antibodies with compound specificities indicated a biochemical relation between P and other blood group systems.

The currently recognized antigens of the P System include two common antigens, P₁ and P, and P_k which has been considered a rare erythrocyte antigen. Approximately 75% of white individuals possess the P₁ antigen, and the P antigen is present on all normal red cells, the erythroblasts and endothelium cells⁷.

Many blood group antigens have been found in numerous parasites, nevertheless the clinical significance of these antigenic similitudes between a parasite and its host is not very clear established yet¹².

The antigens of the P System have also been detected in various animal species: *Taenia Echinococcus*, *Ascaris suum*, *Fasciola gigantica*, *Paragonimus westermani*, *Dicrocoelium dentriticum*¹³. Another excellent source for the P_k antigen is the fluid in hydatid cysts¹.

The blood group P substances of tissue cells mediate the adhesion

of a series of pathogenic bacterias (*Escherichia coli*, *Streptococcus suis*, *Pseudomona aeruginosa* and the Shiga toxin of *Shigella dysenteriae*) and P globoside has been shown to be the cellular receptor for parvovirus B₁₉⁹. At the present time it is accepted that many microorganims and parasites carry blood group antigens which have become evident by Inhibition Agglutination Tests and Immune Fluorescence Techniques. Yet, the clinical significance of this fact is still unknown¹³.

The aim was to study the presence of P₁ antigenic determiners in *Ascaris lumbricoides* extracts by means of the use of different monoclonal antibodies of well-known concentrations and Ig class.

We worked with 14 *A. lumbricoides* extracts [AE]. In order to perform our experiments [AE] were prepared. Adult specimens were washed in physiological solution supplemented with 200 mg/ml of streptomycin and 200 mg/ml of penicillin. After that a refrigerated mechanical rupture was performed for 5 days. The supernatants were collected and kept at -20 °C with a final concentration of timerozal 1:1000^{8,9}. Inhibition Agglutination Test was made facing the [AE] against anti P₁ monoclonal antibodies in optimal concentrations. Suspensions of fresh red cells (P₁) were used as a revealing system⁶. The experiments were carried out in a bromelin enzymatic medium (hydrolysis C-terminal peptide bond of Lys, Ala, Tyr and Gly) because the P System antigens are resistant to the action of this enzyme, their reaction capacity being increased. The temperature chosen was 4 °C because the optimal thermic is room temperature or lower¹¹. The anti P₁ monoclonal antibodies were supplied

by Monoclonal Antibodies against Blood Group Antigens (Monoclonal IV Workshop, Paris, July, 2001), and the Ig concentration in the antibodies was studied by enzyme-linked immunoassays⁵. Titre, Score (S), and Sensitivity Parameter (SP) were determined for each anti P₁ monoclonal antibody against revealing system (P₁ Red cells)^{2,10,15}.

Anti P₁ monoclonal antibodies characteristics are shown on Table 1.

Table 1

Characteristics of the anti P₁ monoclonal antibodies used in the Inhibition Agglutination Test

Anti P ₁ antibody	Ig Class	Concentration (µg/ml)	Titre	S	SP
2.88	Ig M	< 0.5	1/ 256	45	1.02
2.89	Ig M	< 0.5	1/ 256	27	0.63
2.90	Ig A	5.51	1/ 64	24	0.23
2.91	Ig M	87.2	1/ 1024	43	2.78
2.92	Ig M	< 0.5	1/ 256	39	0.85
2.93	Ig M	< 0.5	1/ 16	13	0.06

Ten [AE] inhibited the agglutination of all anti P₁ monoclonal antibodies with red cells (P₁) showing the epitopes present in the [AE] reacted against the six monoclonal antibodies studied.

The 4 remaining [AE] evidenced a different behaviour against anti P₁ monoclonal antibodies. Only some of them were inhibited by the extracts. Inhibition Agglutination Test results of these 4 [AE] are shown on Table 2.

Table 2

Results of the Inhibition Agglutination Test for the extracts with no conventional behavior

[AE] /ANTI P ₁	2.88	2.89	2.90	2.91	2.92	2.93
[AE] 1	NI	NI	I	I	I	I
[AE] 2	NI	NI	I	I	I	I
[AE] 3	NI	I	I	I	I	I
[AE] 4	NI	NI	I	NI	I	I

I = Inhibition Agglutination; NI = No Inhibition Agglutination.

It is demonstrated that all the [AE] have P₁ activity, nevertheless it can only be made evident by means of a monoclonal antibody battery. The [AE] activity is independent of the Ig class, antibody Ig concentration, titre, Score and Sensitivity Parameter. The activity depends on the epitope at which the monoclonal antibody is directed. The use of these antibodies for routine or investigation works leads to standardized methods and a better definition of the specificity of antigen-antibody.

RESUMEN

Presencia de determinantes antigénicas P₁ en *Ascaris lumbricoides*

Los antígenos del Sistema P han sido detectados en numerosos parásitos, bacterias y virus, aunque todavía se desconoce su significado

clínico. El objetivo fue estudiar la presencia de determinantes antigénicas P₁ en extractos de *A. lumbricoides* mediante el uso de 6 anticuerpos monoclonales de concentraciones y clase de Ig conocidas. Se trabajó con 14 extractos de *A. lumbricoides*. Se realizó la prueba de Inhibición de la Aglutinación en medio enzimático de bromelina y temperatura de 4 °C. Se determinó el título, Score y Parámetro de Sensibilidad de cada anticuerpo monoclonal frente a la suspensión de glóbulos rojos usada como sistema revelador. Diez extractos inhibieron la aglutinación de todos los anticuerpos monoclonales anti P₁ con la suspensión de glóbulos rojos P₁. Los 4 extractos restantes sólo inhibieron la aglutinación de algunos de ellos. Se demuestra que todos los extractos estudiados tienen actividad P₁. Esta actividad es independiente del título, Score, Parámetro de Sensibilidad, concentración y clase de Ig del anticuerpo monoclonal y depende del epítopo al cual esta dirigido el anticuerpo.

REFERENCES

- CAMERON, G.L. & STAVELEY, J.M. - Blood group P substance in hydatid cyst fluid. *Nature (Lond.)*, 179: 147-148, 1957.
- GOUEMAND, M. & DELMAS-MARSALET, Y. - *Elements d'Immuno-Hématologie*. Paris, Flammarion, 1967.
- JOHNSTONE, A.P. & TURNER, M.W. - *Immunochemistry 1*. Oxford, Oxford University Press, 1997.
- LANDSTEINER, K. & LEVINE, P. - Further observations on individual differences on human blood. *Proc. Soc. exp. Biol.*, 24: 941-942, 1927.
- LEVINE, P.; BOBBITT, O.B. & WALLER, R.K. - Isoimmunization of a new blood factor in tumor cells. *Proc. Soc. exp. Biol.*, 77: 403-405, 1951.
- MARCELLI, A.; FINE, J.M. & RIBAT, I. - *Techniques in Immunohematologie*. Paris, Flammarion, 1981.
- MARCUS, D.; KUNDU, S. & SUZUKI, A. - The P blood group system: recent progress in Immunochimistry and Genetics. *Sem. Hemat.*, 18: 63-71, 1981.
- MONROY OSTRIA, A.; GOMEZ GUTIERREZ, L.J.; RAMIREZ RAMIREZ, A. & CARRILLO LANDIN, G. - Reconocimiento por inmunotransferencia de antígenos de *Taenia solium* y su larva. *Rev. Lat.-amer. Microbiol.*, 34: 33-38, 1992.
- PONCE DE LEÓN, P. & VALVERDE, J. - P system epitopes in *Ascaris lumbricoides*. *Rev. Inst. Med. trop. S. Paulo*, 44: 115-116, 2002.
- PONCE DE LEÓN, P.; VALVERDE, J. & ZDERO, M. - Aloantibodies ABO in patients with ascariasis. *Rev. Inst. Med. trop. S. Paulo*, 42: 297-298, 2000.
- REID, M. & LOMAR FRANCIS, C. - *The blood group antigen*. London, Academic Press, 1997.
- SALMON, Ch. - *Les groupes sanguins ou l'écriture des genes*. Paris, Masson, 1997.
- SALMON, Ch.; CARTRON, J.P. & ROUGER, P.H. - *Les groupes sanguins chez l'homme*. Paris, Masson, 1991.
- SANGER, R. - An association between the P and Jay systems of blood groups. *Nature (Lond.)*, 176: 1163-1164, 1955.
- VALVERDE DE RASSIA, J. & RASSIA, J. - Soluciones de baja fuerza iónica en las pruebas de Coombs para la detección de anticuerpos anti Rh. *Acta bioquím. clin. lat-amer.*, 16: 295-305, 1982.

Received: 5 November 2002

Accepted: 18 November 2002